FORNO DA LABORATORIO TTN

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*Cruciale per TTN16
Crucible for TTN16
Ø 130 mm*

LABORATORY KILN TTN model

The kiln is built with a steel structure painted on fire at 180 ° C.

The thermal insulation consists of high-density ceramic fiber panels, suitable for the working temperature of the oven.

Heating is obtained with molybdenum disilicide elements, which do not require a protective atmosphere and do not give rise to aging processes.

This has the advantage of a simple and safe operation of the oven and a long life of the resistances with the possibility of replacing even just one of them without altering the behavior of the others.

It is a fixed crucible laboratory kiln suitable for melting a quantity of mixture which, depending on the model chosen, can vary from +/- 600/1000 gr. up to 3000/5000 gr.

The mixture to be melted is placed in the crucible inside the oven, through a hole of suitable size in the upper part of the oven.

As the mixture is melted through the crucible and from the kiln, first, through a suitably sized hole in the bottom of the crucible and then a larger hole in the bottom of the furnace, under the crucible, which runs through the entire insulation of the oven.

COMAND PANEL

The temperature and firing cycle are controlled by a Lumel RE 82 microprocessor programmer.


With this type of programmer it is possible to configure and memorize a maximum of 15 programs each consisting of a maximum of 15 ramps.

TTN-16-5S

*Cruciale per TTN-16-5S
Crucible for TTN-16-5S
Ø 220 mm*

TECHNICAL CHARACTERISTICS

Mod.	Vol [Lt]	Temp. max	Internal dimensions [mm]			External dimensions[mm]			Power kW	V + N	Weight [kG]
			Width [b]	Depth [a]	Height [c]	Width [B]	Depth [A]	Height [C]			
TTN-16	6	1600 °C	170	170	200	850	650	1620	5	230	230
TTN-16-S5	10,5	1600 °C	Ø 260		200	1250	1350	1710	20	400	476