PRESSA DA LABORATORIO MIGNON S - MIGNON SS

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LABORATORY PRESS MIGNON S - MIGNON SS

Built in a compact and elegant structure of steel painted with epoxy paints cooked at 180°C, they are fully hydraulic and are suitable for the production of samples in the laboratory.

The pressing force is exerted through a hydraulic cylinder placed on the upper part of the press and on which the appropriate punch is installed. The mold, placed on the work surface, is manually filled with the powder to be pressed and after pressing the sample obtained is extracted by a piston placed in the lower part.

The pressing force is manually adjustable and is displayed by means of a pressure gauge placed in the upper front part of the machine.

CONTROL PANEL

The front panel, to operate with the machine, presents:

- indicator light (shows if the machine is on or off)
- main switch
- two-handed control system
- extraction button

Supplied with:

- screed
- allen keys for mold assembly / disassembly


Not included:

- mold, supplied at the customer's choice

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
Molde available on request:
special molde Ø 40 mm.
with sintered Zirconium Ox punch

NOTE
The moulds currently supplied are suitable for
use with a maximum pressure on the shape of
1700 Kg/cm²

IMPORTANT
For proper press operation, the pressure range of the circuit at
which you can operate must be within between:
- 30 bar and 120 bar (MIGNON-S)
- 30 bar and 210 bar (MIGNON-SS)

Calculation example:
Piston diameter = 150 mm
Piston area= 176,63 cm²
Shape area = 50 cm²
Desired pressure on the shape = 400 Kg/cm²
Circuit pressure: (400 x 50)/176,63 = 113,23 bar
(pressure value to be set)