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## RAPID MILL WITH PLANETARY DRIVE, MOD. MAG/3000/D WITH 4 ROTATING UNITSS

The requirement of some laboratories in the ceramic industry to be able to rapidly grind large frit quantities, led Ceramic Instruments to design an innovative system that has achieved truly interesting results.

The "MAG" series laboratory mill, uses a planetary-type rotating system where 4 grinding containers (jars) are placed on a bearing plate which rotates in the opposite direction. The energy stored due to the centrifugal forces involved, drags the grinding bodies (alumina balls) and obtain an exceptional grinding effect.

### The main advantages of this system are:

- Up to 6-8 Kg of dry product in one single grinding operation.
- Grinding residue better than that obtained from traditional systems with eccentric rotation.
- Shorter grinding times (40-45 minutes for 6 Kg of ceramic frit with a residue of 1 g).
- Grinding evenness (practically no residual coarse parts).
- Vibration-free and perfect machine stability.
- Less noisy than traditional systems (approx. 72-74 dB).
- Overall dimensions are approx. 1/4 of those of traditional systems with the same production capacity (equal to at least to 3 SD series mills with 4 rotating units each using 1000 cc jars. In figures, 0.99 m2 instead of 4.23 m2).
- Energy saving (it works with one 1.5 kW motor instead of at least 12 0.22 kW motors required for the traditional systems with the same production capacity).

### General technical features:

- Epoxy painted steel sheet bearing structure.
- Height adjustable feet.
- Nicasil treated jar-holding units, other drive components galvanised or burnished.
- Safety micro-switch on door.
- Door lock so that it cannot be opened before the rotating units completely stop.
- Sound-proofing obtained by special sound-absorbing material, placed on the basket protecting the rotating units, on the door and under the drive gears.
- Rotating units that can easily be moved manually when the machine is stopped in order to obtain the ideal position for jar loading and unloading.

### Electrical control panel composed of:

- Main switch 0/1
- Operator panel with 4" touch screen display

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#### Main parameters set by the user:

- Grinding time
- Speed of rotation of containers
- Possibility to reverse rotation
- Reversal rotation time
- Interval time between reversals of rotation
- Possibility of saving 20 customized working programs

Speed of rotation of jars: from 60 to 260 rpm  
Power input: 230 V, single-phases, 50/60 Hz, 1.5 KW.  
Overall dimensions: 113x80x110 cm.  
Net weight 300 Kg

#### Supplied complete with:

four 3 l porcelain jars and relative alumina balls.  
Charge of dry product for each jar 1.5÷2 Kg.  
Total charge 6÷8 Kg.

#### Important:

The MAG/3000/D mills must always work with the units perfectly balanced.  
The jars must always be 4 and must have the same weight.

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## ACCESSORIES AND SPARE PARTS:

### Accessories for emptying jars:

- Stainless steel screen,  $\varnothing$  300 mm height 90 mm.
- Mesh with reinforcing support, net opening 1.00 mm, ASTM 18. Net weight 1.25 kg.
- Stainless steel screen,  $\varnothing$  300 mm, h. 90 mm. Mesh with 0.120 mm net opening, ASTM 120. Net weight 1.1 kg.  
As above, net opening on request.
- Stainless steel connector,  $\varnothing$  300 mm, height 155 mm. Capacity 11 l. Net weight 1.6 kg.

### HARD PORCELAIN JARS:

Made from alumina porcelain fired at a high temperature, totally non-porous. The external surface is glazed white. Mohs Hardness 8.

Specific weight: 2.4 g/cm<sup>3</sup>. AL<sub>2</sub>O<sub>3</sub>: content 40%.

- Supplied complete with cap and seal.
- **3 liters jar.**  
Dimensions:  $\varnothing$  187 mm, 272 mm high.  
Net weight 8,15 kg.
- Charge of alumina balls for above. Total weight 3 kg. Composition: 2,5 kg of  $\varnothing$  18.5 mm balls.  
0.5 kg of  $\varnothing$  25.0 mm balls.
- Seal for 3l jar in rubber
- Knob for cup of 3l jar
- Complete Cap for 3l jar
- **1.5 liters jar.**  
Dimensions:  $\varnothing$  150 mm, 260 mm high  
Net weight 4,5 kg.
- Charge of alumina balls for above. Total weight 1.5 kg.  
Composition: 1 kg of  $\varnothing$  18.5 mm balls.  
0.5 kg of  $\varnothing$  25.0 mm balls.
- Seal 1,5l jar in rubber.
- Knob for cup of 1,5l jar
- Complete Cap for 1,5l jar

### ALUMINA BALLS (Al<sub>2</sub>O<sub>3</sub>):

- Al<sub>2</sub>O<sub>3</sub> content 90%.
  - Specific weight 3.57 g/cm<sup>3</sup>.
  - Mohs hardness 9
  - Porosity: 0.
  - Colour: white.
- 1 kg of balls,  $\varnothing$  18.5 mm. 1 kg of balls,  $\varnothing$  25.0 mm.
  - 1 kg of balls,  $\varnothing$  25.0 mm.
  - 1 kg of balls (microspheres)  
with mixed diameters of 6 mm and 7 mm.

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